

# IL2 Antibody, FITC conjugated



PACO27607

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Produced by T-cells in response to antigenic or mitogenic stimulation, this protein is required for T-cell proliferation and other activities crucial to regulation of the immune response. Can stimulate B-cells, monocytes, lymphokine-activated killer cells, natural killer cells, and glioma cells.

**Gene ID:**

IL2

**Uniprot**

P60568

**Synonyms:**

Interleukin-2 (IL-2) (T-cell growth factor) (TCGF) (Aldesleukin), IL2

**Immunogen:**

Recombinant Human Interleukin-2 protein (21-153AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A