

# IL-33 Antibody



PACO27617

---

## Product Information

**Size:**

50ug

**Reactivity:**

Pig, Mouse, Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000

**Protein Background:**

**Gene ID:**

IL-33

**Uniprot**

M5B263

**Synonyms:**

**Immunogen:**


Recombinant Pig Interleukin-33 protein (123-276AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



### Western blot

All lanes: IL-33 antibody at 2 $\mu$ g/ml

Lane 1: Hela whole cell lysate

Lane 2: Mouse lung tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 32 kDa

Observed band size: 32 kDa