

IL-33 Antibody, Biotin conjugated



PACO27620

Product Information

Size:

50ug

Reactivity:

Pig

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

IL-33

Uniprot

M5B263

Synonyms:

Immunogen:

Recombinant Pig Interleukin-33 protein (123-276AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A