

IL4 Antibody, Biotin conjugated



PACO27624

Product Information

Size:

50ug

Reactivity:

Dog

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Participates in at least several B-cell activation processes as well as of other cell types. It is a costimulator of DNA-synthesis. It induces the expression of class II MHC molecules on resting B-cells. It enhances both secretion and cell surface expression of IgE and IgG1. It also regulates the expression of the low affinity Fc receptor for IgE (CD23) on both lymphocytes and monocytes.

Gene ID:

IL4

Uniprot

O77762

Synonyms:

Interleukin-4 (IL-4) (B-cell stimulatory factor 1) (BSF-1) (Lymphocyte stimulatory factor 1), IL4

Immunogen:

Recombinant Dog Interleukin-4 protein (25-132AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A