

ILVBL Antibody, Biotin conjugated



PACO27644

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

membrane.

Gene ID:

ILVBL

Uniprot

A1L0T0

Synonyms:

Acetolactate synthase-like protein (EC 2.2.1. -) (IlvB-like protein), ILVBL, AHAS

Immunogen:

Recombinant Human Acetolactate synthase-like protein (45-344AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A