

# IPO13 Antibody, FITC conjugated



PACO27659

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

IPO13

**Uniprot**

O94829

**Synonyms:**

Importin-13 (Imp13) (Karyopherin-13) (Kap13) (Ran-binding protein 13) (RanBP13), IPO13, KIAA0724 RANBP13

**Immunogen:**

Recombinant Human Importin-13 protein (501-700AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A