

PACO27661

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcription factor involved in the induction of interferons IFNA and INFB and inflammatory cytokines upon virus infection. Activated by TLR7 or TLR8 signaling.

**Gene ID:**

Irf5

**Uniprot**

P56477

**Synonyms:**

Interferon regulatory factor 5 (IRF-5), Irf5

**Immunogen:**

Recombinant Mouse Interferon regulatory factor 5 protein (1-497AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A