

Itgal Antibody, FITC conjugated



PACO27686

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Integrin alpha-L/beta-2 is a receptor for ICAM1, ICAM2, ICAM3 and ICAM4. Is involved in a variety of immune phenomena including leukocyte-endothelial cell interaction, cytotoxic T-cell mediated killing, and antibody dependent killing by granulocytes and monocytes. Mice expressing a null mutation of the alpha-L subunit gene demonstrate impaired tumor rejection and impaired leukocytes recruitment.

Gene ID:

Itgal

Uniprot

P24063

Synonyms:

Integrin alpha-L (CD11 antigen-like family member A) (Leukocyte adhesion glycoprotein LFA-1 alpha chain) (LFA-1A) (Leukocyte function-associated molecule 1 alpha chain) (Lymphocyte antigen 15) (Ly-15) (CD antigen CD11a), Itgal, Lfa-1 Ly-15

Immunogen:

Recombinant Mouse Integrin alpha-L protein (153-325AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A