

Itgb1 Antibody, Biotin conjugated



PACO27695

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Itgb1

Uniprot

P09055

Synonyms:

Integrin beta-1 (Fibronectin receptor subunit beta) (VLA-4 subunit beta) (CD antigen CD29), Itgb1

Immunogen:

Recombinant Mouse Integrin beta-1 protein (21-728AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A