

KIAA0513 Antibody



PACO27781

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:2000-1:5000,
IF:1:50-1:200

Protein Background:**Gene ID:**

KIAA0513

Uniprot

O60268

Synonyms:

Uncharacterized protein KIAA0513, KIAA0513


Immunogen:

Recombinant Human Uncharacterized protein KIAA0513 protein (1-301AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot

Positive WB detected in: Jurkat whole cell lysate, K562 whole cell lysate, Mouse brain tissue

All lanes: KIAA0513 antibody at 2.7 μ g/ml


Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 47, 34, 46 kDa

Observed band size: 57 kDa

.



Immunofluorescent analysis of HepG2 cells using PACO27781 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).