

KIAA0513 Antibody, Biotin conjugated



PACO27784

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	KIAA0513
Source:	Uniprot
Rabbit	O60268
Isotype:	Synonyms:
IgG	Uncharacterized protein KIAA0513, KIAA0513
Applications:	Immunogen:
ELISA	Recombinant Human Uncharacterized protein KIAA0513 protein (1-301AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A