

Klk7 Antibody, Biotin conjugated



PACO27820

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Glandular kallikreins cleave Met-Lys and Arg-Ser bonds in kininogen to release Lys-bradykinin. Predominant kallikrein protein in the kidney.

Gene ID:

Klk7

Uniprot

P36373

Synonyms:

Glandular kallikrein-7, submandibular/renal (rGK-7) (EC 3.4.21.35) (Esterase B) (Kallikrein-related protein K1) (Proteinase A) (RSKG-7) (Tissue kallikrein), Klk7, Klk-7

Immunogen:

Recombinant Rat Glandular kallikrein-7, submandibular/renal protein (25-261AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A