

Kng1 Antibody, FITC conjugated



PACO27839

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

Kng1

Uniprot

P08934

Synonyms:

Kininogen-1 [Cleaved into: Kininogen-1 heavy chain; Bradykinin; Kininogen-1 light chain], Kng1, Kng

Immunogen:

Recombinant Rat Kininogen-1 protein (2-110AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A