

# Kng1 Antibody, Biotin conjugated



PACO27840

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Kng1

**Uniprot**

P08934

**Synonyms:**

Kininogen-1 [Cleaved into: Kininogen-1 heavy chain; Bradykinin; Kininogen-1 light chain], Kng1, Kng

**Immunogen:**

Recombinant Rat Kininogen-1 protein (2-110AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A