

KRT17 Antibody



PACO27853

Product Information

Size:

50ug

Reactivity:

Bovine

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

KRT17

Uniprot

A1L595

Synonyms:

Keratin, type I cytoskeletal 17 (Cytokeratin-17) (CK-17) (Keratin-17) (K17), KRT17

Immunogen:

Recombinant Bovine Keratin, type I cytoskeletal 17 protein (1-441AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A