

# KRT17 Antibody, Biotin conjugated



PACO27856

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Bovine	KRT17
<b>Source:</b>	<b>Uniprot</b>
Rabbit	A1L595
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Keratin, type I cytoskeletal 17 (Cytokeratin-17) (CK-17) (Keratin-17) (K17), KRT17
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Bovine Keratin, type I cytoskeletal 17 protein (1-441AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A