

Lgals3bp Antibody, FITC conjugated



PACO27939

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Lgals3bp

Uniprot

O70513

Synonyms:

Galectin-3-binding protein (Cyp-C-associated protein) (CyCAP) (Lectin galactoside-binding soluble 3-binding protein) (Mac-2-binding protein) (MAC2BP) (Mac-2 BP) (Protein 90K), Lgals3bp

Immunogen:

Recombinant Rat Galectin-3-binding protein (19-574AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A