

# Lhcgr Antibody, FITC conjugated



PACO27946

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Receptor for lutropin-choriogonadotropic hormone. The activity of this receptor is mediated by G proteins which activate adenylate cyclase.

**Gene ID:**

Lhcgr

**Uniprot**

P16235

**Synonyms:**

Lutropin-choriogonadotropic hormone receptor (LH/CG-R) (Luteinizing hormone receptor) (LSH-R), Lhcgr

**Immunogen:**

Recombinant Rat Lutropin-choriogonadotropic hormone receptor protein (27-362AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A