

Lhcgr Antibody, Biotin conjugated



PACO27948

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor for lutropin-choriogonadotropic hormone. The activity of this receptor is mediated by G proteins which activate adenylate cyclase.

Gene ID:

Lhcgr

Uniprot

P16235

Synonyms:

Lutropin-choriogonadotropic hormone receptor (LH/CG-R) (Luteinizing hormone receptor) (LSH-R), Lhcgr

Immunogen:

Recombinant Rat Lutropin-choriogonadotropic hormone receptor protein (27-362AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A