

LIG1 Antibody, Biotin conjugated



PACO27960

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

DNA ligase that seals nicks in double-stranded DNA during DNA replication, DNA recombination and DNA repair.

Gene ID:

LIG1

Uniprot

P18858

Synonyms:

DNA ligase 1 (EC 6.5.1.1) (DNA ligase I) (Polydeoxyribonucleotide synthase [ATP] 1), LIG1

Immunogen:

Recombinant Human DNA ligase 1 protein (1-919AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A