

# Lilrb3 Antibody, HRP conjugated



PACO27966

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May act as receptor for class I MHC antigens. Becomes activated upon coligation of LILRB3 and immune receptors, such as FCGR2B and the B-cell receptor. Down-regulates antigen-induced B-cell activation by recruiting phosphatases to its immunoreceptor tyrosine-based inhibitor motifs (ITIM).

**Gene ID:**

Lilrb3

**Uniprot**

P97484

**Synonyms:**

Leukocyte immunoglobulin-like receptor subfamily B member 3 (LIR-3) (Leukocyte immunoglobulin-like receptor 3) (Cell-surface glycoprotein p91) (Paired immunoglobulin-like receptor B) (PIR-B), Lilrb3, Pirb

**Immunogen:**

Recombinant Mouse Leukocyte immunoglobulin-like receptor subfamily B member 3 protein (664-841AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A