

Lilrb3 Antibody, Biotin conjugated



PACO27968

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May act as receptor for class I MHC antigens. Becomes activated upon coligation of LILRB3 and immune receptors, such as FCGR2B and the B-cell receptor. Down-regulates antigen-induced B-cell activation by recruiting phosphatases to its immunoreceptor tyrosine-based inhibitor motifs (ITIM).

Gene ID:

Lilrb3

Uniprot

P97484

Synonyms:

Leukocyte immunoglobulin-like receptor subfamily B member 3 (LIR-3) (Leukocyte immunoglobulin-like receptor 3) (Cell-surface glycoprotein p91) (Paired immunoglobulin-like receptor B) (PIR-B), Lilrb3, Pirb

Immunogen:

Recombinant Mouse Leukocyte immunoglobulin-like receptor subfamily B member 3 protein (664-841AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A