

LOXL1 Antibody



PACO27989

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:10-1:100

Protein Background:

Active on elastin and collagen substrates.

Gene ID:

LOXL1

Uniprot

Q08397

Synonyms:

Lysyl oxidase homolog 1 (EC 1.4.3. -) (Lysyl oxidase-like protein 1) (LOL), LOXL1, LOXL


Immunogen:

Recombinant Human Lysyl oxidase homolog 1 protein (106-574AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO27989 at dilution of 1:20.