

# LOXL1 Antibody, HRP conjugated



PACO27991

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Active on elastin and collagen substrates.

**Gene ID:**

LOXL1

**Uniprot**

Q08397

**Synonyms:**

Lysyl oxidase homolog 1 (EC 1.4.3. -) (Lysyl oxidase-like protein 1) (LOL), LOXL1, LOXL

**Immunogen:**

Recombinant Human Lysyl oxidase homolog 1 protein (106-574AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A