

LRRTM2 Antibody



PACO28029

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the development and maintenance of excitatory synapse in the vertebrate nervous system. Regulates surface expression of AMPA receptors and instructs the development of functional glutamate release sites. Acts as a ligand for the presynaptic receptors NRXN1-A and NRXN1-B.

Gene ID:

LRRTM2

Uniprot

O43300

Synonyms:

Leucine-rich repeat transmembrane neuronal protein 2 (Leucine-rich repeat neuronal 2 protein), LRRTM2, KIAA0416 LRRN2

Immunogen:

Recombinant Human Leucine-rich repeat transmembrane neuronal protein 2 protein (165-422AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A