

LY6H Antibody



PACO28062

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

nervous system development, organ morphogenesis.

Gene ID:

LY6H

Uniprot

O94772

Synonyms:

Lymphocyte antigen 6H (Ly-6H), LY6H


Immunogen:

Recombinant Human Lymphocyte antigen 6H protein (26-115AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry analysis of human testis tissue using PACO28062 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO28062 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry analysis of human tonsil tissue using PACO28062 at dilution of 1:100.