

# Lypla1 Antibody, HRP conjugated



PACO28091

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Hydrolyzes fatty acid, from S-acylated cysteine residues in proteins such as trimeric G alpha proteins or HRAS. Has depalmitoylating activity toward KCNMA1. Has low lysophospholipase activity.

**Gene ID:**

Lypla1

**Uniprot**

P97823

**Synonyms:**

Acyl-protein thioesterase 1 (APT-1) (EC 3.1.2. -) (Lysophospholipase 1) (Lysophospholipase I) (LPL-I) (LysoPLA I), Lypla1, Apt1 Pla1a

**Immunogen:**

Recombinant Mouse Acyl-protein thioesterase 1 protein (1-230AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A