

MAP2K6 Antibody, FITC conjugated



PACO28144

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MAP2K6

Uniprot

P52564

Synonyms:

Dual specificity mitogen-activated protein kinase kinase 6 (MAP kinase kinase 6) (MAPKK 6) (EC 2.7.12.2) (MAPK/ERK kinase 6) (MEK 6) (Stress-activated protein kinase kinase 3) (SAPK kinase 3) (SAPKK-3) (SAPKK3), MAP2K6, MEK6 MKK6 PRKMK6 SKK3

Immunogen:

Recombinant Human Dual specificity mitogen-activated protein kinase kinase 6 protein (1-334AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A