

# MAS1L Antibody, HRP conjugated



PACO28183

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytoplasm, integral component of membrane, integral component of plasma membrane, nucleoplasm, plasma membrane, G-protein coupled receptor activity, G-protein coupled receptor signaling pathway.

**Gene ID:**

MAS1L

**Uniprot**

P35410

**Synonyms:**

Mas-related G-protein coupled receptor MRG (MAS-R) (MAS1-like), MAS1L, MRG

**Immunogen:**

Recombinant Human Mas-related G-protein coupled receptor MRG protein (318-378AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A