

MAS1L Antibody, Biotin conjugated



PACO28185

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoplasm, integral component of membrane, integral component of plasma membrane, nucleoplasm, plasma membrane, G-protein coupled receptor activity, G-protein coupled receptor signaling pathway.

Gene ID:

MAS1L

Uniprot

P35410

Synonyms:

Mas-related G-protein coupled receptor MRG (MAS-R) (MAS1-like), MAS1L, MRG

Immunogen:

Recombinant Human Mas-related G-protein coupled receptor MRG protein (318-378AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A