

MAT2A Antibody, Biotin conjugated



PACO28193

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the formation of S-adenosylmethionine from methionine and ATP.

Gene ID:

MAT2A

Uniprot

P31153

Synonyms:

S-adenosylmethionine synthase isoform type-2 (AdoMet synthase 2) (EC 2.5.1.6) (Methionine adenosyltransferase 2) (MAT 2) (Methionine adenosyltransferase II) (MAT-II), MAT2A, AMS2 MATA2

Immunogen:

Recombinant Human S-adenosylmethionine synthase isoform type-2 protein (1-395AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A