

Mcl1 Antibody, FITC conjugated



PACO28200

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in the regulation of apoptosis versus cell survival, and in the maintenance of viability but not of proliferation. Mediates its effects by interactions with a number of other regulators of apoptosis. Isoform 2 has antiapoptotic activity.

Gene ID:

Mcl1

Uniprot

P97287

Synonyms:

Induced myeloid leukemia cell differentiation protein Mcl-1 homolog (Bcl-2-related protein EAT/mcl1), Mcl1

Immunogen:

Recombinant Mouse Induced myeloid leukemia cell differentiation protein Mcl-1 homolog protein (1-308AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A