

MFSD11 Antibody



PACO28230

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

MFSD11

Uniprot

O43934

Synonyms:

UNC93-like protein MFSD11 (Major facilitator superfamily domain-containing protein 11) (Protein ET), MFSD11, ET


Immunogen:

Recombinant Human UNC93-like protein MFSD11 protein (191-238AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO28230 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung tissue using PACO28230 at dilution of 1:100.