

# MIF4GD Antibody, HRP conjugated



PACO28243

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Functions in replication-dependent translation of histone mRNAs which differ from other eukaryotic mRNAs in that they do not end with a poly-A tail but a stem-loop. May participate in circularizing those mRNAs specifically enhancing their translation.

**Gene ID:**

MIF4GD

**Uniprot**

A9UHW6

**Synonyms:**

MIF4G domain-containing protein (SLBP-interacting protein 1) (hSLIP1), MIF4GD, SLIP1

**Immunogen:**

Recombinant Human MIF4G domain-containing protein (1-133AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A