

# MLLT10 Antibody, FITC conjugated



PACO28284

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Probably involved in transcriptional regulation. In vitro or as fusion protein with KMT2A/MLL1 has transactivation activity. Binds to cruciform DNA.

**Gene ID:**

MLLT10

**Uniprot**

P55197

**Synonyms:**

Protein AF-10 (ALL1-fused gene from chromosome 10 protein), MLLT10, AF10

**Immunogen:**

Recombinant Human Protein AF-10 protein (1-126AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A