

Mmp7 Antibody



PACO28298

Product Information

Size:

50ug

Reactivity:

Rat, Mouse, Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:**Gene ID:**

Mmp7

Uniprot

P50280

Synonyms:

Matrilysin (EC 3.4.24.23) (Matrin) (Matrix metalloproteinase-7) (MMP-7) (Pump-1 protease) (Uterine metalloproteinase), Mmp7, Mmp-7


Immunogen:

Recombinant Rat Matrilysin protein (98-267AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: Mmp7 antibody at 2 μ g/ml

Lane 1: Colo320 whole cell lysate

Lane 2: NIH/3T3 whole cell lysate

Lane 3: Mouse lung tissue

Lane 4: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

.