

# MRPS5 Antibody



PACO28386

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000

**Protein Background:**

mitochondrial inner membrane, mitochondrial small ribosomal subunit, mitochondrion, poly(A) RNA binding, structural constituent of ribosome, mitochondrial translation, mitochondrial translational elongation, mitochondrial translational initiation, mitochondrial translational termination, organelle organization.

**Gene ID:**

MRPS5

**Uniprot**

P82675

**Synonyms:**

28S ribosomal protein S5, mitochondrial (MRP-S5) (S5mt) (Mitochondrial small ribosomal subunit protein uS5m), MRPS5

**Immunogen:**

Recombinant Human 28S ribosomal protein S5, mitochondrial protein (161-430AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



### Western blot

All lanes: MRPS5 antibody at 5 $\mu$ g/ml

Lane 1: HeLa whole cell lysate

Lane 2: Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 49, 28 kDa

Observed band size: 49 kDa

.