

# CKB Antibody, Biotin conjugated



PACO28393

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Reversibly catalyzes the transfer of phosphate between ATP and various phosphogens (e. g. creatine phosphate). Creatine kinase isoenzymes play a central role in energy transduction in tissues with large, fluctuating energy demands, such as skeletal muscle, heart, brain and spermatozoa.

**Gene ID:**

CKB

**Uniprot**

P12277

**Synonyms:**

Creatine kinase B-type (EC 2.7.3.2) (Brain creatine kinase) (B-CK) (Creatine kinase B chain) (Creatine phosphokinase B-type) (CPK-B), CKB, CKBB

**Immunogen:**

Recombinant Human Creatine kinase B-type protein (2-270AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A