

PACO28406

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Membrane-anchored forms may play a role in cellular adhesion. Megakaryocyte-potentiating factor (MPF) may potentiate megakaryocyte colony formation.

Gene ID:

Msln

Uniprot

Q61468

Synonyms:

Mesothelin (Pre-pro-megakaryocyte-potentiating factor) [Cleaved into: Megakaryocyte-potentiating factor (MPF); Mesothelin, cleaved form], Msln, Mes Mpf

Immunogen:

Recombinant Mouse Mesothelin protein (298-600AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A