

# MTHFR Antibody, FITC conjugated



PACO28412

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the conversion of 5,10-methylenetetrahydrofolate to 5-methyltetrahydrofolate, a co-substrate for homocysteine remethylation to methionine.

**Gene ID:**

MTHFR

**Uniprot**

P42898

**Synonyms:**

Methylenetetrahydrofolate reductase (EC 1.5.1.20), MTHFR

**Immunogen:**

Recombinant Human Methylenetetrahydrofolate reductase protein (1-656AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A