

MTHFR Antibody, Biotin conjugated



PACO28417

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the conversion of 5,10-methylenetetrahydrofolate to 5-methyltetrahydrofolate, a co-substrate for homocysteine remethylation to methionine.

Gene ID:

MTHFR

Uniprot

P42898

Synonyms:

Methylenetetrahydrofolate reductase (EC 1.5.1.20), MTHFR

Immunogen:

Recombinant Human Methylenetetrahydrofolate reductase protein (1-656AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A