

MTMR11 Antibody



PACO28426

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable pseudophosphatase. Contains a Glu residue instead of a conserved Cys residue in the dsPTPase catalytic loop which renders it catalytically inactive as a phosphatase (Potential).

Gene ID:

MTMR11

Uniprot

A4FU01

Synonyms:

Myotubularin-related protein 11 (Cisplatin resistance-associated protein) (hCRA), MTMR11

Immunogen:

Recombinant Human Myotubularin-related protein 11 protein (1-240AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A