

PACO28438

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Prothymosin alpha may mediate immune function by conferring resistance to certain opportunistic infections.

**Gene ID:**

PTMA

**Uniprot**

P06454

**Synonyms:**

Prothymosin alpha [Cleaved into: Prothymosin alpha, N-terminally processed; Thymosin alpha-1], PTMA, TMSA

**Immunogen:**


Recombinant Human Prothymosin &alpha; protein (5-110AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO28438 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO28438 at dilution of 1:100.