

# MXD3 Antibody, Biotin conjugated



PACO28445

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transcriptional repressor. Binds with MAX to form a sequence-specific DNA-binding protein complex which recognizes the core sequence 5'-CAC[GA]TG-3'. Antagonizes MYC transcriptional activity by competing for MAX and suppresses MYC dependent cell transformation.

**Gene ID:**

MXD3

**Uniprot**

Q9BW11

**Synonyms:**

Max dimerization protein 3 (Max dimerizer 3) (Class C basic helix-loop-helix protein 13) (bHLHc13) (Max-associated protein 3) (Max-interacting transcriptional repressor MAD3) (Myx), MXD3, BHLHC13 MAD3

**Immunogen:**

Recombinant Human Max dimerization protein 3 protein (1-206AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A