

# Myc Antibody



PACO28454

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Myc

**Uniprot**

P01108

**Synonyms:**

Myc proto-oncogene protein (Proto-oncogene c-Myc) (Transcription factor p64), Myc

**Immunogen:**

Recombinant Mouse Myc proto-oncogene protein (1-439AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A