

NAPRT Antibody, Biotin conjugated



PACO28497

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the conversion of nicotinic acid, (NA) to NA mononucleotide (NaMN).
Essential for NA to increase cellular NAD levels and prevent oxidative stress of the cells.

Gene ID:

NAPRT

Uniprot

Q6XQN6

Synonyms:

Nicotinate phosphoribosyltransferase (NAPRTase) (EC 6.3.4.21) (FHA-HIT-interacting protein) (Nicotinate phosphoribosyltransferase domain-containing protein 1), NAPRT, FHIP NAPRT1

Immunogen:

Recombinant Human Nicotinate phosphoribosyltransferase protein (1-310AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A