

# NARS Antibody, Biotin conjugated



PACO28501

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytoplasm, cytosol, extracellular exosome, asparagine-tRNA ligase activity, gene expression, tRNA aminoacylation for protein translation.

**Gene ID:**

NARS

**Uniprot**

O43776

**Synonyms:**

Asparagine--tRNA ligase, cytoplasmic (EC 6.1.1.22) (AsparaginyI-tRNA synthetase) (AsnRS), NARS

**Immunogen:**

Recombinant Human Asparagine--tRNA ligase, cytoplasmic protein (1-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A