

NDUFA6 Antibody



PACO28522

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

NDUFA6

Uniprot

P56556

Synonyms:

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 6 (Complex I-B14) (CI-B14) (LYR motif-containing protein 6) (NADH-ubiquinone oxidoreductase B14 subunit), NDUFA6, LYRM6 NADHB14


Immunogen:

Recombinant Human NADH dehydrogenase [ubiquinone] 1 α subcomplex subunit 6 protein (27-154AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO28522 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung cancer using PACO28522 at dilution of 1:100.