

NDUFA6 Antibody, Biotin conjugated



PACO28525

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NDUFA6

Uniprot

P56556

Synonyms:

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 6 (Complex I-B14) (CI-B14) (LYR motif-containing protein 6) (NADH-ubiquinone oxidoreductase B14 subunit), NDUFA6, LYRM6 NADHB14

Immunogen:

Recombinant Human NADH dehydrogenase [ubiquinone] 1 α subcomplex subunit 6 protein (27-154AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A