

PACO28526

Product Information

Size:

50ug

Reactivity:

Staphylococcus aureus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Carrier of the growing fatty acid, chain in fatty acid, biosynthesis. Is able to confer high methicillin resistance to S.aureus when overproduced.

Gene ID:

acpP

Uniprot

P0A001

Synonyms:

Acyl carrier protein (ACP), acpP, hmrB

Immunogen:

Recombinant Staphylococcus aureus Acyl carrier protein (1-77AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A